



Micro-USB Connector



To avoid safety hazards and short circuits, DO NOT connect both batteries into the PCB at the same time.

LIPO Battery Charger




U1 is the battery charger that charges either the auxillary LiPo battery (labelled VBAT\_JST2) or the coin cell Li-Ion battery connected directly into Biocoin (labelled VBAT\_JST1).

To charge the auxillary LiPo battery, short pins 1 and 2 on the J2 jumper (labelled VBAT\_AUX on PCB).

To charge the coin cell Li-Ion battery, short pins 2 and 3 on the J2 jumper (labelled VBAT\_COIN on PCB).


M2 and J3 allow VBAT\_JST2 to power Biocoin if USB power is not connected, or allows USB +5V to power Biocoin if connected and also disconnects VBAT\_JST2 from Biocoin to allow the battery to be charged without simultaenously discharging.



Title		
Size	Number	Revision
B		
Date:	10/02/2025	Sheet of
File:	C:\Users\...\POWER.SchDoc	Drawn By:


Digital/Power POGO Pin Signals

These two sets of pins contain the digital signals and power nets that connect up to Biocoin for debug, programming, charging, etc.




SWD Header

Digital SWD signals needed for burning the MCU bootloader or programming the device via Segger




Digital Signal Headers

Headers for digital signal debugging/observation




Title		
Size	Number	Revision
B		
Date:	10/02/2025	Sheet of
File:	C:\Users\...\DIGITAL_DATA.SchDoc	Drawn By:

# Analog POGO Pin Signals




J11 and J13 pins contain the analog signals coming from external electrode connectors and send them up to Biocoin.




# Electrode Connectors

SPE Connector 1




Supports either a 3-electrode connector or a 4-electrode connector. Uses WE1 by default and also WE2 if using 4-electrode connector.

SPE Connector 2



3 electrode measurement for 1 of the OCP measurement channels.

# Electrode Connectors




Phoenix screw terminal connectors. Sends all 13 electrode signals to these pins + GND to enable connection to other types of sensors or electrode cables.



# Dummy Loads


Mock cells for mimicking different types of electrochemical/impedance measurements.



Dummy loads for BioZ:  
- Component values picked to mimic those in the AN-1557 App Note from Analog Devices

For 2-wire BioZ:  
- Use CE0 and AIN1, and disconnect jumpers to AIN2 and AIN3

For 4-wire BioZ:  
- Use all 4 signals and populate all jumpers



Left load: Dummy load for voltammetry/amperometry. To easily test chronoamperometry, don't populate 1M resistor

Right load: Dummy load for EIS. To easily test EIS, don't populate the capacitor

In software, must switch the MUX to WE4 to use dummy loads here

Title		
Size	Number	Revision
B		
Date:	10/02/2025	Sheet of
File:	C:\Users\...\ANALOG_DATA.SchDoc	Drawn By: